



**DISCLAIMER** :This document was produced to support MSF staff taking into account MSF implementation conditions and practices and is intended for internal use only. This document is based on MSF's experience and current knowledge which makes it provisional, not necessarily applicable to all contexts nor useful for all users, and subject to regular revisions and changes, if not removal altogether. MSF Operational Centre Brussels does not accept any responsibility or liability with regard to the use of information contained in this document by any third party.





## List of necessary equipment for oxygen therapy

---


*Annex 4 - Oxygen therapy and Oxygen humidification*



## List of necessary equipment for oxygen therapy

### Annex 4 - Oxygen therapy and Oxygen humidification

| Code        | Label  | Notes/Pictures  |
|-------------|--|---|
| SCTDMASO1A- | OXYGEN FACE MASK, simple, with tubing, adult size            |    |
| SCTDMASO1P- | OXYGEN FACE MASK, simple, with tubing, pediatric size        |   |
| SCTDMASOH1A | OXYGEN MASK, HIGH CONCENTRATION, adult size                  |    |
| SCTDMASOH1P | OXYGEN MASK, HIGH CONCENTRATION, pediatric size              |   |
| SCTDCANN2PL | NASAL OXYGEN CANNULA, 2 prongs + tube, premature low flow    | Premature low-flow nasal cannulas are for neonates below 2.5kgs.<br> |
| SCTDCANN2N- | NASAL OXYGEN CANNULA, 2 prongs + tube, neonate               |   |
| SCTDCANN2P- | NASAL OXYGEN CANNULA, 2 prongs + tube, pediatric             |   |
| SCTDCANN2A- | NASAL OXYGEN CANNULA, 2 prongs + tube, adult                 |   |
| SCTDCONN7A1 | CONNECTOR, biconical, asym., ext Ø 7-11 & 10-13mm, autoclav  |    |
| SCTDCONN7S- | CONNECTOR, biconical, symmetrical, ext. Ø 7-11 mm, autoclav. |   |
| SCTDTUBE052 | TUBE, silicone, autoclavable, int. Ø 5 mm, 25 m              | Can be cut to shorter lengths to avoid tripping and coiling of the tube.  |
| EEMDCONE4-- | CONCENTRATOR O2 (DeVilbiss 525KS) 5l, 220V + access.         | 5 litre/min max   |

|             |  |                  |
|-------------|--|------------------|
| EEMDCONE5-- | CONCENTRATOR O2 (New Life Intensity) 10l, 230V, 50 Hz + acc. | 10 litre/min max |
|-------------|--|------------------|

|             |  |   |
|-------------|--|---|
| EEMDCONE6-- | CONCENTRATOR O2 portable (Eclipse 5), 3l + acc.              | 3 litre/min max   |
| EEMDCONS335 | (conc. DeVilbiss 515KS/AKS/525KS) OUTLET CONNECTOR plast.gre |    |
| EEMDCONA003 | HUMIDIFIER, autoclavable (FLOWMETER TR200)                   |    |
| EEMDCONE010 | KIT FLOW SPLITTER pediatric (Sureflow) + HUMIDIFIERS         |    |
| EEMDCONA008 | (ped. distr. Sureflow) EXTRA INTERCONNECTION HOSE, HO151-1   |  |
| EEMDCONA009 | (ped. distr. Sureflow) INTERCONNECTION HOSE, 9/16 F 10meters |   |
| EEMDCONA010 | (ped. distr. Sureflow) INTERCONNECTION HOSE, 9/16 F 25meters |   |
| EEMDCONA011 | (ped. distr. Sureflow) INTERCONNECTION HOSE, 9/16 F 15meters |   |
| EEMDCONA012 | (ped. distr. Sureflow) INTERCONNECTION HOSE, 9/16 F 20meters |   |
| EEMDCONA008 | (ped. distr. Sureflow) EXTRA INTERCONNECTION HOSE, HO151-1   |   |

|               |  |   |
|---------------|--|---|
| SDREEAHC3S- - | DRESSING HYDROCOLLOID, thin<br>for neonate / malnourished patient / sensitive skin |  |
| SDRETAPA1F25  | TAPE, ADHESIVE, ELASTIC FOAM, roll, 2.5 cm x 5 m                                   |   |
| DEXTALCO5S-   | ALCOHOL-BASED HAND RUB, solution, 500 ml, bot.                                     |   |
| EEMDCONC313   | (conc. DeVilbiss 515KS/AKS/525KS) ANTIBACT. FILTER PV5LD-651                       |   |
| KMEDMTUB1--   | (oxygen concentrator 5-10 litres) MODULE TUBERCULOSIS                              |  |
| SDISSUQA2B-   | DETERGENT/DISINFECTANT for surfaces, 2 l tin + dosing pump                         | Or the equivalent solution available within the project                             |