



DISCLAIMER :This document was produced to support MSF staff taking into account MSF implementation conditions and practices and is intended for internal use only. This document is based on MSF's experience and current knowledge which makes it provisional, not necessarily applicable to all contexts nor useful for all users, and subject to regular revisions and changes, if not removal altogether. MSF Operational Centre Brussels does not accept any responsibility or liability with regard to the use of information contained in this document by any third party.



When to use nasal cannula, a simple mask or a non-rebreathing mask?

Annex 1 - Oxygen therapy and Oxygen humidification

When to use nasal cannula, a simple mask or a non-rebreathing mask?

Annex 1 - Oxygen therapy and Oxygen humidification

Type:	Nasal Cannula/Prongs
When to use:	When delivery low-flow oxygen (0.1L-4L in paediatrics and up to 6L in adults) Preferred method of oxygen delivery in children less than 5 years of age.
Advantages:	Come in varying sizes from preterm to adult (4 in total) More comfortable than masks and helps with compliance In paediatrics: Higher concentrations may be given via nasal cannula
Disadvantages:	Higher flow rates can irritate and dry the nasal mucosa. If high-flow oxygen therapy via nasal cannula to be given over 2 hours, humidified oxygen should be used In neonates and infants, improperly sized cannula can cause nasal obstruction or epistaxis Nasal cannulas are contraindicated in patients with a nasal deformity or trauma, surgery, epistaxis, fractured base of skull, nasal/sinus congestion Frequent clearing of the nasal passages may be required to optimize efficacy of nasal cannula (see procedure on nasopharyngeal suctioning if necessary) CAUTION: Nasal Cannula/Prongs have an increased risk of aerosolization compared to simple oxygen masks or NRM.
Type:	Simple Mask
When to use:	When delivering oxygen above 5L/min between in paediatrics and between 5L/min-10L/min in adults
Advantages:	Can be used for patients experiencing nasal irritation or epistaxis Useful for patients who are strictly mouth breathers
Disadvantages:	Available in only 2 sizes (paediatric and adult) and may not fit well to each face Less comfortable than nasal cannula, muffles communication and can cause increased distress/anxiety in patients. Restricts drinking and feeding CAUTION: If used with a lower flow rate, the risk of inhaling CO ₂ that has been exhaled into the mask increases, thus, increases the risk of altered conscious states (hypercapnia) and worsen respiratory distress. CAUTION: Ensure each mask is tightly fitted, to reduce risk of aerosol spread
Type:	Non-Rebreathing/high concentration
When to use:	For flow rates over 5L/min in paediatrics and over 8L/min in patients over 12 years old* *The healthcare provider must be capable of ensuring an adequate and maintained flow sufficient to ensure the reservoir bag remains 2/3rds full during inspiration. Ideally, the oxygen flow should be above a flow of oxygen ≥ 10 -12 L/min. Currently, the MSF oxygen concentrators can only provide a maximum flow rate of 10L/min. When possible, flow rates should be greater than the above values.
Advantages:	80-100% oxygen can be administered Suitable for patients in respiratory distress or any type of shock
Disadvantages:	Available in only 2 sizes (paediatric and adult) and may not fit well to each face The healthcare provider must ensure a good seal between the mask and the patient and Less comfortable than nasal cannula, muffles communication and can cause increased distress/anxiety in patients. Restricts drinking and feeding

	 <p>CAUTION: If the oxygen flow rate is under 5L/min in paediatrics or under 8L/min in adults, the bag could be partly filled by CO₂ exhaled by the patient. Rebreathing his/her own CO₂ could alter conscious state (hypercapnia) and worsen the respiratory distress</p>
	 <p>CAUTION: Ensure each mask is tightly fitted, to reduce risk of aerosol spread, and do not remove the little green valves from the NRM.</p>